




e•SPEC™ Products

	PerkinElmer	Spectron (std)	e•SPEC™ NEW!
<u>Pt Sampler</u>	WE027802	SC2013-Pt \$2,050	SC2013-Pt_e \$1,800
<u>Pt Skimmer</u>	WE027803	SC2014-Pt \$1,450	SC2014-Pt_e \$1,200



e•SPEC platinum cones for the Elan 6000 and 9000 series ICP-MS


	Agilent	Spectron (std)	e•SPEC™ NEW!
<u>Pt Sampler</u> <i>18mm Pt tip</i>	G1820-65360	HP1006A-Pt \$3,600	HP1006A-Pt_e \$3,000
<u>Pt Sampler</u> <i>10mm Pt tip</i>	G1820-65239	HP1006-Pt \$1,650	HP1006-Pt_e* \$1,300
<u>Pt Skimmer</u>	G1833-65132	HP1008CS-Pt \$1,350	HP1008CS-Pt_e \$1,175



e•SPEC HP1006-Pt_e sampler
Displayed with Pt and Ni skimmers,
Pt shields for 4500 and 7500

*Standard 10mm Pt tip in a 25mm nickel core for greater longevity over the Pt in Cu design.


	Thermo Finnigan	Spectron (std)	e•SPEC™ NEW!
<u>Sampler</u>	1044530	T1001-Ni \$225	T1001TE-Ni_e \$250
<u>Skimmer</u>	1067600	T1002A-Ni \$225	T1002TE-Ni_e \$250
<u>Pt Sampler</u> <i>15.2 mm Pt tip</i>	1067500 <i>6mm tip</i>	T1006-Pt/Cu \$1,750	T1006TE-Pt_e \$1,500
<u>Guard Electrode</u>	1226640	**	T5004 \$700



e•SPEC Thermally Enhanced (TE)
cones for Element 2 & Neptune in nickel
(standard) and aluminum (optional), with
Pt Guard Electrode

High TDS cones for the Thermo Finnigan Element 2 and Neptune series ICP-MS

	Thermo Elemental Xi Cones	Spectron (std) Xi Cones	e•SPEC™ NEW!
<u>Sampler</u>	3600812	VG1021-Ni/Cu \$290	VG1021TE-Ni/Cu_e \$ 290
<u>Skimmer</u>	3600811	VG1022-Ni \$200	**



e•SPEC Thermally Enhanced (TE)
cones for X-Series ICP-MS systems

Spectron Xi cones are available in nickel/copper, nickel and aluminum configurations

*Delivery outside the United States will be influenced by local taxes